

PRODUCT MOEDL



500/500 LED
Adjust color temperature



700/500 LED
Adjust color temperature



700/700 LED
Adjust color temperature



500 LED
Adjust color temperature



700 LED
Adjust color temperature

700/500 SERIES TECHNICAL DATA

Type	Illuminance	Colour temperature	Colour reduction index (Ra)	Total irradiance	Size of light field	Service life of illuminant	Power supply voltage	Best height of installation	Total power consumption	Total LED bulb quantity	ENDO
500 / 500 LED	≥140,000Lux ≥140,000Lux	3800±500K 4400±500K 5000±500K	93	534 W/m ² 534 W/m ²	250-300mm	50,000h	AC110-240V 50/60Hz	2900mm	60W 60W	144pcs (12*6+12*6)	YES
700 / 500 LED	≥140,000Lux ≥160,000Lux	3800±500K 4400±500K 5000±500K	93	544 W/m ² 534 W/m ²	250-300mm	50,000h	AC110-240V 50/60Hz	2900mm	60W 80W	180pcs (18*6+12*6)	YES
700 / 700 LED	≥160,000Lux ≥160,000Lux	3800±500K 4400±500K 5000±500K	93	544W/m ² 544W/m ²	250-300mm	50,000h	AC110-240V 50/60Hz	2900mm	80W 80W	218pcs (18*6+18*6)	YES
500 LED	≥140,000Lux	3800±500K 4400±500K 5000±500K	93	534W/m ²	250-300mm	50,000h	AC110-240V 50/60Hz	2900mm	60W	72pcs (12*6)	YES
700 LED	≥160,000Lux	3800±500K 4400±500K 5000±500K	93	544W/m ²	250-300mm	50,000h	AC110-240V 50/60Hz	2900mm	80W	108pcs (18*6)	YES

Two color opts for the mask in the lamp head

