

Features

- The device adopts design of foreign advanced model advantages, using micro-computer control, high performance and reliable quality with automatic constant temperature, washing and film output.
- The device is wall mounted with two types of common and heating. The heating one has special function of automatic constant temperature, which ensures the reliability and high quality when developing the film with the poor conditions in cold weather.
- Light room operation, small volume, save space, low noise, operation easy and fast. Replace the cumbersome process of traditional hand-washing film, which overcome problems of low speed and instable quality.
- Developing various oral films in 2-8minutes, and it can be adjusted according to the exposure time.
- ☐ Can use various films size2x3、3x4 or 2.7x3.4cm.
- Indication function of liquid expire date, the device is equipped with memory disk of liquid replaced date, convenient the customer to replace liquid timely.





Autamatic heats installment

Dark Box

- It is suitable for the development of all sizes of dental films and occlusive films(not include panoramic films). It instead of darkroom because of the shield against light. The box is light weight, convenient placement and carrying.
- There are 3 cups in the box,put the developer,fixer and clear water in turn ,the liquid height should over the height of the films.

