

Anesthesia System & Ventilator Specialist

Main Features

- ◆ Working mode: Continuous
- ◆ Medical Air compressor
- ♦ Working pressure display
- ◆ Compact design, easy to transfer
- Brake of front wheels
- ◆ Thomas pump, quality assurance
- ◆ Can be equipped with international ventilator brand



Main Specifications

Output pressure: 0.4MPa±0.1MPa

output flow ≥45L/min

Peak flow: When duration is 1 second, the peak flow is

≥120L/min

Dew point temperature: ≤ 5°C

Alarm temperature: ≥123°C±10°C

Assembled with clean control module, supplying cleaner and

safer air to the patients

Power: 650VA

Noise: ≤50dB



