

CE

EXCEED IN QUALITY AND MORE

BA. 680
Urine Analyzer



BA. 680

Urine Analyzer

General Features

Throughput: 500 samples/hour
Principle: Dual wavelength reflectance photometry
Touch screen
Sensor for auto detecting urine strip
Easy maintenance design with auto waste collection
Support auto&manual print
Original urine strips



Technical Specifications

Parameters: Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leukocyte, Vitamin C
Input: 5.6-inch touch screen
Output: Internal thermal printer
Display: LCD
Lamp source: LED, 470nm, 525nm, 625nm
Storage: Up to 2,000 results
Working Temperature: 10-30°C
Humidity: 20%-80%
Dimension: 390*330*210mm
Power: AC100~240V, 50/60 Hz, 45W
Net Weight :6 kg

Easy operation

User-friendly software allows you a better experience on the analyzer
Intelligent sensor make urine strip detected automatically

Less maintenance

Auto waste collection and alert make it less contaminated as well as less maintenance
Power cut off protection make it even less break-down

High efficiency

High throughput with 500 samples/hour, ideal for extensive sample laboratory
Auto print result after the test is finished

