

HF A-1 NATURAL COLOR ELASTIC CREPE BANDAGE

Material: 80% cotton; 20% spandex

Weight gsm: 75g,80g,90g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm, 7.5cm, 10cm, 15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.





HF A-4 RED(BLUE) THREAD ELASTIC CREPE BANDAGE

Material: 80% cotton; 20% spandex Weight gsm: 70g,80g,100g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer's Request.



HF A-2 SKIN COLOR ELASTIC CREPE BANDAGE

Material: 80% cotton; 20% spandex Weight gsm: 62g,70g,75g,80g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm, 7.5cm, 10cm, 15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



**HF A-5 COTTON CREPE BANDAGE (100% COTTON)** 

Material: 100% pure cotton Weight gsm: 75g,80g,85g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane Note: Personalized specifications as possible at

customer's Request.



HF A-3 BLEACHED ELASTIC CREPE BANDAGE

Material: 80% cotton; 20% spandex Weight gsm: 62g,70g,75g,80g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm, 7.5cm, 10cm, 15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



### HF B-1 B-3 NATURAL(BLEACHED) COLOR PLAIN ANDAGE

Material: 55% pure cotton;45% spandex

Weight gsm: 75g,80g,85g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m

Width: 5cm, 7.5cm, 10cm, 15cm

Basic packing: Individually packed in cellophane Note: Personalized specifications as possible at

customer's Request.



HF B-2 SKIN COLOR ELASTIC PLAIN BANDAGE

Material: 55% cotton; 45% spandex

Weight gsm: 75g,80g,90g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



HF D-1 GAUZE BANDAGE

Material: 100% cotton Yarn: 40's, 32's, 21's

Mesh: 19x15,30x20,30x30,30x28,40x20,28x24 Size:Length(stretched):4m,4.5m,5m,6m,9m,10m Width: 4cm,5cm,6cm,7.5cm,8cm,10cm,15cm Basic packing: Individually packed in cellophane Note: Personalized specifications as possible at

customer's Request.



HF C-1 BLEACHED TABBY PBT BANDAGE

Material: 55% viscose; 45% cotton with woven edges

Weight gsm: 28g,30g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



HF F-1 POP PLASTER OF PARIS BANDAGE

Material: cotton and plaster of paris

Color: white

Solidifying time: 2 to 5 minutes

Size:Length(stretched):2.7m,3m,3.6m,4m,4.6m,5m Width: 5cm,7.5cm,10cm,12.5cm,15cm,20cm

Basic packing: Individually packed in waterproof paper

Note: Personalized specifications as possible at

customer's Request.



HF C-6 THICK PBT BANDAGE

Material: 55% viscose; 45% cotton with woven edges

Weight gsm: 28g,30g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



#### HF F-2 ORTHOPEDIC BANDAGE

Material: cotton covered nonwoven,PP wool

Color: white

Weight gsm: 90g,135g

Size:Length(stretched):2.7m,3m,4m

Width: 5cm,7.5cm,10cm,12.5cm,15cm,20cm Basic packing: Individually packed in cellophane Note: Personalized specifications as possible at

customer's Request.



HF F-3 ELASTIC COHESIVE BANDAGE

Material: Elastic fiber and cotton fiber Weight gsm: 30g,34g,78g,80g

Color: Vaious

Size:Length(stretched):4m,4.5m,5m

Width: 2.5cm,5cm,7.5cm,10cm,15cm,20cm Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



HF B2 HIGH ELASTIC BANDAGE

Material:75% polyester;25% artificial silk

Weight gsm: 100g,120g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m

Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane Note: Personalized specifications as possible at

customer's Request.



HF A1 SKIN COLOR HIGH ELASTIC BANDAGE

Material:70% Skin color polyester;30% rubber

Weight gsm: 90g,100g,110g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



HF B2 BLUE THREAD HIGH ELASTIC BANDAGE

Material:50% cotton;20% polyester,30% rubber

Weight gsm: 90g,100g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at

customer's Request.



HF B1 HIGH ELASTIC BANDAGE

Material:85% silver pink polyester;15% rubber

Weight gsm: 81g,90g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



#### HF F1 SKIN COLOR HIGH ELASTIC BANDAGE

Material:57% cotton;26.5% polyester,16.5% rubber

Weight gsm: 80g,85g,90g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m

Width: 5cm,7.5cm,10cm,15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at

customer's Request.



HF H1 HIGH ELASTIC BANDAGE

Material:50% Cotton;20% polyester;30% rubber

Weight gsm: 100g,110g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm, 7.5cm, 10cm, 15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer

's Request.



HF K1 PURFLED HIGH ELASTIC BANDAGE

Material:55% Cotton;20% polyester;25% rubber

Weight gsm: 90g,100g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm, 7.5cm, 10cm, 15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer Basic packing: Individually packed in cellophane

's Request.



### HF J-1 FIRST AID BANDAGE BANDAGE

Material:Bandage and Pad

Bandage: Conforming bandage or gauze bandage

Pad: Low adherent dressing, cotton covered

with gauze or nonwoven Size:Bandgage LxW(Pad LxW)

6cmx4m(6x8cm);8cmx4m(8x10cm)

10cmx4m(10x12cm);12cmx5m(12x16cm)

Note: Personalized specifications as possible at customer's Request.



HF N1 PURFLED HIGH ELASTIC BANDAGE

Material:52% Cotton;20% spandex;38% fiber

Weight gsm: 80g,90g,100g

Clip: With or without clips, elastic band clips, metal clips

Size:Length(stretched):4m,4.5m,5m Width: 5cm, 7.5cm, 10cm, 15cm

Basic packing: Individually packed in cellophane

Note: Personalized specifications as possible at customer Length: 10cmx1.75cm(Original Length)

's Request.



HF J-2 FIRST AID BANDAGE(H-TYPE)

200pcs/carton; Carton size: 76x55x38cm

Cotton Pad Size: 20x24cm

EO sterilization











HF J-3 GREEN FIRST AID BANDAGE

Length: 10cmx4m(Stretched Length) Cotton Pad Size: 10cmx10cm Sealed Package;200pcs/Carton Carton size: 50x43x43cm

**EO** Sterilization







HF JT01 Cervical Collar:

Shaped for confort and to give the correct support to help relieve the pain suffered from whiplash, stiff neck and muscular tension.

Made from high density closed cell foam(4D&100D)

for better support, have a woven covering. Available in 3size: small, medium and large S:44x9cm M:48x10cm L:53x11cm



HF OS01 Orthopedic Splints(Brace):

Lower leg, Night Spline, Walker Style,

below Knee, Hook and Loop Strapes. Material: nylon, polyester, neoprene.

Size: Various Sizes Color: Back



HF OS02 Orthopedic Splints:

For Arm, Hook and Loop Strapes Material: Nylon, Polyester, Neoprene Size: Various Color: Black











Cotton Stretch Bandage/Crepe: made by pure cotton;108gsm

Elasticity: The fully-streched length is not less than twice the unstretched length; the regin length is not more than Two-Thirds the fully-stretched length.

Threads per 10cn Warp: singles 55 to 65, two fold 55 to 65.

Weight per unit area: Not less than 105gsm

Size: 5cmx4.5m;17.5cmx4.5m; 10cmx4.5m: 15cmx4.5m



HF CS01
Casting Tapes:
230gsm;Made by polyester or glass fiber
Size:2"/3"/4"/5"/6" x 4YDS

Packed by box



HF W01 Wasit Support, Magnetic back support

Materil: Nylon,Polyester,Neoprene and Magnet Size: S,M,L,XL Color: Black



HF Z-1
Cuff, made by pure cotton and foam;
Size: 5cmx6m; 7,5cmx6m
Package:2/1 rolls into one box



HF Z-2 Sterotrac --- Skin Traction Kit

Size	Ref	Pack
Adhesive Adult	9515	1x18
Adhesive Child	9516	1x24
Non-Adhesive Adult	9517	1x18
Non-Adhesive Child	9518	1x24



HF Z-3 Vaseline Gauze

10cmx10cm,

1pc/Sterile Package; 10packs/box

Gamma Sterilization

# Vaseline Bandage

10cmx10m;10cmx8m

1roll/Sterile Package; Gamma Sterilization



HF Z-4 Tubular bandage

Besides the above specifications, we can also produce other specifications according the samples and the customer's requirement.

Specification: 2"/3"/4" x 10m/25m



# HF Z-5

## Zinc oxide tape

- 1. Sawtooth zinc oxide plaster available
- 2. Different size and packing available
- 3. Main material: natural rubber and zinc oxide, 100% cotton
- 4. Scope of Application:surgical&sport using

Specification: 2"/3"/4" x 5m/5yards/10m/10yards



## HF Z-6 Triangular bandage

- 1, Material:,100% cotton or Nonwoven fabric
- 2, Certificate:CE, FDA,ISO
- 3, Different size is available

Specification: 90x90x127cm, 96x96x136cm 100x100x142cm