





Brief Introduction

■ Effective, User-friendly, Intelligent transmission

3.0 inch color touch screen, large character indicates and read measurement value clearly and rapidly;

On-screen display of multiple measurement data, and simultaneous storage of multiple patient information, achieving data management to patients.

■ Long endurance, environment friendly and durable

Equipped with charging base, enduring built-in chargeable battery, measuring 2000 times after battery is full;

Xenon flash lamp with long lifespan (measurement ≥150,000 times);

Charging base with the function of integrated light checker, consumable-free, saving cost for check.

Accurate measurement and quick check

Effective check to identify patient with high-risk jaundice precisely, conform to the newest clinical guidance;

Constant measurement to multiple patients without repeating calibration; Average measurement: this function permits 1~5 times setting of average measurement to ensure the accuracy on results;

Converting between two kinds of measurement unit mg/dL, µmol/L.

To the second state of the

Integrated light checker

Standard Configuration

Consist of main body, base and power adapter.

