




# Cytotoxic Safety Cabinet



**LCD Display**  
Large display is easy to monitor all the safety parameters at a glance, and ergonomically sized control panel improves users interface.




**Water and Gas Tap**  
High quality stainless steel




**Front Window**  
Motorized, Two-layer laminated toughened glass >5mm, anti UV



**Foot Switch**  
Adjust front window height by foot during experiment, to avoid airflow turbulence caused by arm movement.



**UV Lamp**  
Emission of 253.7 nanometers for most efficient decontamination



**Bag In and Bag Out Filter**  
The PVC bag can ensure the safety of the personnel and the environment, and make the replacement process convenient and quick.

## TECHNICAL PARAMETERS

Model	11234BBC86
Internal Size (W*D*H)	1220*600*630mm 48" *23.6"*24.8"
External Size (W*D*H)	1370 *760 *2100 mm 53.9"*29.9" *86.6"
Tested Opening	Safety Height ≤ 200mm
Max Opening	400mm
Internal Work Area, space	0.73M <sup>2</sup>
Inflow Velocity	0.46m/s
Down Flow Velocity	0.33m/s
Exhaust Volume with Thimble Duct	CBV Exhaust Volume 611 m <sup>3</sup> / h (360 cfm) Static Pressure at CBV Exhaust Volume 39 Pa / 0.15 in H <sub>2</sub> O
Air Supply Filter	ULPA filter, efficiency > 99.999% for particle size between 0.1 -0.2 um
Exhaust Filter	2 pieces ULPA filter, efficiency > 99.999% for particle size between 0.1 -0.2 um Filter life indicator and Filter working time indicator.
Alarm	Audio and visual alarm, filter replacement, Window over height, Abnormal airflow velocity
Noise	≤ 60dB
Water & Gas Tap	Water tap*1; Gas tap*1 30W*1
UV Lamp	Emission of 253.7 nanometers for most efficient decontamination; UV life indicator
Fluorescent Lamp	28W*2
Illumination	≥ 1000Lux
Consumption	≤ 700W
Waterproof Socket	Two, Total load of two sockets: ≤ 500W
Display	LCD display
Control System	Microprocessor
Humidity and Temperature	Humidity and Temperature can be displayed on the control panel
Material	Work area: 304 stainless steel Main body: cold-rolled steel sheets with anti-bacteria powder coating
Front Window	Motorized, Two-layer laminated toughened glass ≥ 5mm, anti UV
Caster	Universal caster with brake and leveling feet
Power Supply	110/220V ± 10%, 60/50Hz
Standard Accessory	Base stand, Drain valve, waterproof sockets*2, Remote control, Foot switch, UV lamp*2 and Fluorescent Lamp, Water tap, Gas tap
Optional Accessory	Armrest, Formalin fumigation sterilizer, Infrared sterilizer, Airflow tester
Gross Weight	350kg
Package Size	1570*950*2220mm