



CM-1800

Flat Screen Monitor Chart

Fast Accurate Steady

Features

19-inch high contrast and bright display

Smooth and quiet operation


Wireless Remote control

Wide range of test types


Precise distance adjustment

Easy mounting

Random test sequence



Rotating Astigmatic Fan Chart Fast and Accurate test for Cylinder and Axis



CM-1800

Support data connection with CV-7000 computerized vision tester



Specifications

Display	19" LCD (color) High resolution Monitor
Brightness	250 cd/m ²
Chart versions	E, C, ABC, 123, child, special version charts
Chart acuity	6/3-6/150 20/10-20/500 0.04-2.0
Distance	2.0 ~ 7.0m
Power supply	Input: AC110-220V±10%, 50/60Hz
Power consumption	≤45W
Dimension	455mm (L) × 420mm(W) × 110mm(H)
Weight	4.2kg (Net. Weight)
Remote	IR
Batteries for remote control	3A