


>>> EMV-200 Electric upper gastrointestinal endoscopy



- ◆ Use Sony high-resolution CCD with 440,000 pixels color image
- ◆ With 30,000 hours long life LED lam and high color temperature light cold source
- ◆ Germany Elaborated designed operation system
- ◆ Excellent image freezing and playback function
- ◆ With fully Waterproof design

Diameter of distal-end	Φ9.8mm/ 9.2mm
Diameter of insert tube	Φ9.8mm/ 9.2mm
Clamp aperture	Φ2.8mm
Working length	1035mm
Total length	1335mm
View of field	140°
Depth of view	3-100mm
Resolution	440,000 Pixels high resolution Photosensitive CCD
Tip deflection	up 210° down 90° L/R 100°



>>> EMV-400 Electric colonoscopy



Diameter of distal-end	Φ13.0mm
Diameter of insert tube	Φ12.8mm
Clamp aperture	Φ3.2mm/ 3.7mm
Working length	1330mm
Total length	1630mm
View of field	140°
Depth of view	3-100mm
Resolution	440,000 Pixels high resolution Photosensitive CCD
Tip deflection	up180° down 180° L/R 160°

