



# DENSITY BALANCE

Density balance use high precision weighing sensor and testing device, Its have units convention storage and calculating function.The detection instrument is widely used in scientific research institutions and enterprises

MODEL	XY3002CM	XY6002CM	JA203M	JA602M
Capacity (g)	310	610	210	610
Price(USD)				
Sensor	Load Cell		Load Cell	
Resolution (g)	0.01		0.001	0.01
Min weighing (g)	0.04		0.004	0.04
Pan Size	160x70mm		160x70mm	
Operation temp	5-25°C		5-25°C	
Repeat ability	±0.02g		±0.002g	±0.02g
Liner	±0.03g		±0.003g	±0.03g
Packing Size	440×310×300mm x2		440×310×300mm x2	
Test Mode	Solid、Liquid		Solid、Liquid、Powder、Absorbent etc	
Case Material	Plastic Case with Aluminum density kit		Aluminum Case with Aluminum density kit	
G.weight	4.0kg		6.0 kg	
Display	LCD (white back light with black font)			
Cal.weight	200g			
Tank Size	139x100x78mm			
Option	RS232C & RS485、 Printer			
Unit Mode	Weighing:g/ct/oz	Density: D	Weighing: g/ct/oz	Density: D/%/K...

### Standard Configuration:



### Optional Accessories:



XY-2CM



JA-M