

## YXZ-C501

### Five function electric hospital bed

#### Specifications:

- ◆ Size:L213\*W105\*H48~70cm
- ◆ Back section adjustment range: 0~ 75°±5°
- ◆ Knee section adjustment range: 0~ 40°±5°
- ◆ Height adjustment rang: 20±2cm
- ◆ Trendelenburg & Reverse-Trendelenburg: 0~10°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Epoxy painted steel Structure
- ◆ Wholly molded cold steel sleeping platform,10 pieces divided platform
- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Side rails:Tuck-away side rail,ABS, four pieces full length protection
- ◆ Central locking system with 125mm diameter castor
- ◆ Nurse control panel embedded on the foot board
- ◆ Built-in control embedded panel in the side railing
- ◆ CPR handle accessible from both sides of the bed

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C501B

### Five function electric hospital bed

#### Specifications:

- ◆ Size: L213\*W105\*H48~70cm
- ◆ Backrest section adjustment range: 0~75°±5°
- ◆ Knee section adjustment range: 0~45°±5°
- ◆ Height adjustment range: 20±2cm
- ◆ Trendelenburg & Reverse-Trendelenburg: 0~12°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ One-step punching bed surface. powder coated steel. good air permeability and anti-rust
- ◆ Up-down ABS rails,2pcs/size,safe and reliable
- ◆ Center-control braking system,diameter:150mm

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C502

### Five function electric hospital bed

#### Specifications:

- ◆ Size:L213\*W105\*H48~70cm
- ◆ Back section adjustment range:0~ 75°±5°
- ◆ Knee section adjustment range: 0~40°±5°
- ◆ Height adjustment range:22±2cm
- ◆ Trendelenburg & Reverse-Trendelenburg: 0~12°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Electrophoresis and powder coating
- ◆ Wholly molded cold steel sleeping platform
- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Four pieces of ABS Side rails
- ◆ Central locking system with 125mm diameter castor
- ◆ Platform connected with soft ABS link
- ◆ Built-in control embedded panel in both side rail
- ◆ CPR handle accessible from both sides of the bed, IV pole

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C503

### Five function electric hospital bed

#### Specifications:

- ◆ Size:L213\*W105\*H48~70cm
- ◆ Back adjustment range: 0~80°±5°
- ◆ Knee adjustment range: 0~45°±5°
- ◆ Height adjustment range: 20±2cm
- ◆ Trendelenburg & Reverse-Trendelenburg: 0~12°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Collapsible 4pcs ABS rails,nice appearance
- ◆ Center control braking system,castor diameter:15cm
- ◆ One-step punching bed surface. powder coated steel. good air permeability and anti-rust
- ◆ Equipped with drainage bag hooks.I.V. pole holder, traction pole holder
- ◆ Customized design offered

#### Optional accessories:

- ◆ Movable dinner table, bedside cabinet, Back-up battery, CPR control system, manual crank



# HOSPITAL BED SERIES

## YXZ-C504

### Five function electric hospital bed

#### Specifications:

- ◆ Size:L213\*W105\*H48~70cm
- ◆ Backrest adjustment range:0~75°±5°
- ◆ Knee section adjustment range:0~45°±5°
- ◆ Height adjustable:20±2cm
- ◆ Trendelenburg & Reverse-Trendelenburg:0~12°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Detachable ABS and S.S head and foot board
- ◆ One-step punching bed surface. powder coated steel. good air permeability and anti-rust
- ◆ Up-down ABS rails, 1pc/side, safe and reliable
- ◆ Center-control braking system, diameter:150mm

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C505

### Five function electric hospital bed

#### Specifications:

- ◆ Size:L208\*W96\*H48~70cm
- ◆ Back adjustable:0~80°±5°
- ◆ Knee adjustable:0~45°±5°
- ◆ Height adjustable:23±2cm
- ◆ Trendelenburg & Reverse-Trendelenburg:0~12degree
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Collapsible 2pcs Aluminum alloy rails, nice appearance
- ◆ Rubber covered noiseless castors, all with brake
- ◆ One-step punching bed surface. powder coated steel. good air permeability and anti-rust
- ◆ Equipped with drainage bag hooks, I.V. pole holder, Traction pole holder
- ◆ Customized design offered

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



# HOSPITAL BED SERIES

## YXZ-C506

### Five function electric hospital bed

#### Specification:

- ◆ Size:L213\*W105\*H48~70cm
- ◆ Backrest adjustment range: 0~75°±5°
- ◆ Knee section adjustment range: 0~45°±5°
- ◆ Height adjustment range:20±2cm
- ◆ Trendelenburg & Reverse-Trendelenburg: 0~12°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ ABS bed surface, durable and anti-scratch
- ◆ Up-down ABS rails, 2pc/side, safe and reliable
- ◆ Center-control braking system, diameter:150mm

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C507

### Five function electric hospital sitting bed

#### Specification:

- ◆ Size:L208\*W96\*H48~70cm
- ◆ Backrest adjustment range: 0~80°±5°
- ◆ Knee section adjustment range: 0~45°±5°
- ◆ Height adjustment range:28cm
- ◆ Trendelenburg & Reverse-Trendelenburg: 0~12°
- ◆ Change to chair position: 90°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Collapsible 2pcs Aluminum alloy rails, nice appearance
- ◆ One-step punching bed surface, powder coated steel, good air permeability and anti-rust
- ◆ Rubber covered noiseless castors, diameter 125mm, all with brake
- ◆ Equipped with drainage bag hooks, I.V. pole holder
- ◆ Special design for patient caring

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C301

### Three function electric hospital bed

#### Specifications:

- ◆ Size:L213\*W105\*H48~70cm
- ◆ Backrest section adjustment range:0~75°±5°
- ◆ Knee section adjustment range:0~45°±5°
- ◆ Height adjustment range:20±2cm
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ One-step punching bed surface, durable and anti-scratch
- ◆ Up-down ABS rails, 2pcs/side, safe and reliable
- ◆ Center-control braking system,diameter:150mm

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C302

### Three function electric hospital bed

#### Specifications:

- ◆ Size:L208\*W96\*H48~70cm
- ◆ Backrest section adjustment range:0~75°±5°
- ◆ Knee section adjustment range:0~45°±5°
- ◆ Height adjustment range:20±2cm
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Four section bed platform, powder coating
- ◆ Wholly molded cold steel sleeping platform
- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Two pieces of Aluminum alloy guardrail
- ◆ Remote control panel
- ◆ Luxurious swivel castors: diameter125mm, all with brake
- ◆ Equipped with drainage bag hooks, I.V. pole holder, Traction pole holder
- ◆ Customized design offered

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C303

### Three function electric hospital bed

#### Specifications:

- ◆ Size:L208\*W96\*H48~70cm
- ◆ Backrest section adjustment range:0~75°±5°
- ◆ Knee section adjustment range:0~45°±5°
- ◆ Height adjustment range:20±2cm
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Four section bed platform , powder coating
- ◆ Wholly molded cold steel sleeping platform
- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Two pieces of Aluminum alloy guardrail
- ◆ Remote control panel
- ◆ Center-control braking system: diameter125mm
- ◆ Equipped with drainage bag hooks, I.V. pole holder, Traction pole holder
- ◆ Customized design offered

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C304

### Three function electric hospital bed

#### Specifications:

- ◆ Size:L208\*W93\*H40~70cm
- ◆ Backrest section adjustment range: 0~75°±5°
- ◆ Knee section adjustment range: 0~45°±5°
- ◆ Height adjustment range: 30±2cm
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Four section of mesh bed platform , powder coating
- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Two pieces of Aluminum alloy guardrail
- ◆ Remote control panel
- ◆ Noiseless swivel castors: diameter125mm, all with brake

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



# HOSPITAL BED SERIES

## YXZ-C305

### Three function super-low electric hospital bed

#### Specifications:

- ◆ Size:L208\*W96\*H31~51cm
- ◆ Backrest section adjustment range:0~75°±5°
- ◆ Knee section adjustment range:0~45°±5°
- ◆ Height adjustment range:25cm
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Two pieces of Aluminum alloy guardrail
- ◆ Remote control panel
- ◆ Luxurious swivel castors: diameter125mm, all with brake
- ◆ Equipped with drainage bag hooks, I.V. pole holder, Traction pole holder
- ◆ Customized design offered

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C201

### Two function electric hospital bed

#### Specifications:

- ◆ Size:L208\*W96\*H50cm
- ◆ Backrest section adjustment range:0~75°±5°
- ◆ Knee section adjustment range :0~45±5°±5°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Four section of Mesh bed platform , powder coating
- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ Two pieces of Aluminum alloy guardrail
- ◆ Remote control panel
- ◆ Noiseless swivel castors: diameter125mm, all with brake
- ◆ Equipped with drainage bag hooks, I.V. pole holder, Traction pole holder
- ◆ Customized design offered

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



# HOSPITAL BED SERIES

## YXZ-C202

### Two function electric hospital bed

#### Specifications:

- ◆ Size:L208\*W96\*H50cm
- ◆ Backrest section adjustment range:0~75°±5°
- ◆ Knee section adjustment range :0~45°±5°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Four section of bed platform , powder coating
- ◆ Detachable ABS Head and Foot board, durable and easy to clean
- ◆ One-step punching bed surface, powder coated steel, good air permeability and anti-rust
- ◆ Two pieces of Aluminum alloy guardrail, nice appearance
- ◆ Remote control panel
- ◆ Luxurious swivel castors: diameter125mm, all with brake
- ◆ Equipped with drainage bag hooks, I.V. pole holder, Traction pole holder
- ◆ Customized design offered

#### Optional accessories:

- ◆ Mattress, Movable dinner table, bedside cabinet, Back-up battery, CPR control system



## YXZ-C203

### Electric Two function folding hospital bed

#### Specifications:

- ◆ Size:L200\*W90\*H48.5cm
- ◆ Backrest section adjustment range:0~75°±5°
- ◆ Knee section adjustment range :0~45°±5°
- ◆ Loading capacity: 240kgs

#### Description:

- ◆ Four section strip bed platform ,powder coated steel
- ◆ Easy to fold, portable
- ◆ Durable and powder coated head and foot board, easy to clean.
- ◆ Two pieces of powder coated guardrail, can be dropped
- ◆ Remote control panel
- ◆ Luxurious swivel castors: diameter125mm, all with brake
- ◆ Equipped with drainage bag hooks, I.V. pole holder
- ◆ Customized design offered

#### Optional accessories:

- ◆ Mattress, bedside cabinet, Back-up battery

