

Medical Anti-fog Disposable face Shield

(Isolation face Mask)



Specifications :

Double side Anti-fog, anti-static, No-glare lens;

Light weight and comfortable to wear , quick and easy to don;

Specific Foam band can absorb sweat and provide enough room for eye glasses or safety goggles

Foam design for increase air flow and comfort;

Excellent protection against potential contamination from blood pathogens, body fluids, or harmful chemical splash

Size : 33*22cm

Material : PET+ Sponge +Elastic band

Standard CE EN166, FDA listing

Package : 1pcs/PE Bag

200 Pcs /CTN , 66cm X 35cm X 42cm

G.W.: 8.6kg / N.W.: 7.6kg

CE FDA

