

FFP2 Particulate Protection Face Mask

EN149:2001 + A1:2009, CE 1282 FFP2 NR

Watch its video on YouTube: <https://youtu.be/xrePcw6TU3c>

Company Name:
Company Address:
Tel:
Fax:
Website:

Pioway Medical Lab Equipment Co. Ltd.
Room 1316, De Ying International Plaza, Nanjing, 210012, China
+86-25-84677741
+86-25-84677742
www.pioway.com

Performance:

To prevent the particulate air pollutants, like dust, smoke and smog etc; good protection again viruses and respiratory diseases.



Open the mask, hang the rope on the ear, adjust the position of the bridge of nose, so that the mask covers the mouth and nose and fits the face.

Warning:

1. Failure to follow all instructions and limitations on the use of this product could reduce the effectiveness of respirator and result in sickness or death
2. A properly selected respirator is essential to protect your health
3. This product does not supply oxygen
4. Do not use when concentrations of contaminants are immediately dangerous to health or life
5. Do not use in explosive atmospheres

Parameters:

Product Name: Particulate Respirator

Filter Material: Mechanical and Electrostatic Filter


Protection Level: FFP2



Wearing: Ear wearing type

Nose Clip: Aluminum Bar Clip

Certification: CE Certification, Passed the Latest Revised European Standard EN149: 2001

Pictures:





Certificate:

CERTIFICATE ■ CERTIFICADO ■ CERTIFIKAT



VERIFICATION OF COMPLIANCE

Verification No.:CPAH20031522556

Applicant:
 Address of Applicant: NO1, Hexiang Road, Changfeng, Hefei, Anhui,China
 Manufacturer:
 Address of
 Manufacturer: NO1, Hexiang Road, Changfeng, Hefei, Anhui,China
 Product Name: Disposable mask
 Model No.: M-9501
 Test Standard(s): EN 149: 2001+A1:2009
 Test report(s): TPAH20031522556
 Date of issue: March 15th, 2020
 Date of Expiry: March 14th, 2025

In accordance with the following Applicable Directives:
 The Personal Protective Equipment Regulation (EU)2016/425

Conclusion
 The referred test report(s) show that the product complies with standard(s) recognized as giving presumption of compliance with the essential requirements.
 The CE marking as shown below can be affixed on the product after preparation of necessary technical documentation.
 Other relevant Directives have to be observed.
 Approved by: Hermann Weiher


Global Testing Services Co., Ltd
 E-mail:info@gts-lab.com http://www.gts-lab.com
 Floor 2nd, Building D-1, No. 128, Shehu Road, Minhang District, Shanghai, China.




For and on behalf of
 Global Testing Services Co., Ltd

Copyright of this certificate is owned by Global Testing Services Co., Ltd and may not be reproduced other than for use and with the prior approval of the General Manager.