

H618 HbA1c Meter

Technical specifications

- ▶ Test parameters: HbA1c,eAG
- ▶ Sample type: whole blood
- ▶ Sample volume: 10 μ L
- ▶ Test time: 3-4 min.
- ▶ Measurement range: HbA1c 4.0%-15.0%
- ▶ CV: \leq 5%
- ▶ Storage: 10,000 sample; 12 months QC value, power failure protection
- ▶ COM: RS-232, HIS/LIS
- ▶ Power requirement: AC220V \pm 10%, 50Hz,50VA, with Li battery option
- ▶ Dimension: 205mm \times 130mm \times 70mm, 1KG
- ▶ Working temperature: 5-40 $^{\circ}$ C
- ▶ Consumables: HbA1c cartridge, kits, IC card

CE
ISO13485



H618 HbA1c Meter



- Dry system, maintenance free
- Multi-color touch screen
- Disposable cartridge, point of care