

## FERTILITY



## HCG Pregnancy Test Kits

One Step Cassette Style HCG (human chorionic gonadotropin) Urine Pregnancy Test is a highly sensitive test kit for the determination of HCG in urine (and serum) specimens. This test kit is used to obtain a visual, qualitative result for early detection of pregnancy. After getting the CE certificate, our HCG test kits are widely sold to EU market.

### Features

1. An easy to use pregnancy test. No additional instrument or reagent is needed.
2. Detect within one week of pregnancy;
3. Accuracy: 99.9% Sensitivity: 10mIU/ml, or regular 25 mIU/ml.
4. Designed for both professional and home users.
5. Single package are available for the OTC market.

### Test Procedures

#### Strip



#### Cassette



#### Midstream



### Ordering information

CAT. NO.	Description	Method	Packing	Specimen
HCG 112	HCG Urine Strip	GICA*	100 Tests/Box; 48 Boxes/Carton	Urine
HCG 113	HCG Urine Cassette	GICA*	25 Tests/Box; 40 Boxes/Carton	Urine
HCG 113B	HCG Urine Cassette (Box)	GICA*	25 Tests/Box; 40 Boxes/Carton	Urine
HCG 114	HCG Urine Midstream	GICA*	1 Test/Box; 360 Boxes/Carton	Urine
HCG 115	HCG U/S Combo Strip	GICA*	100 Tests/Box; 48 Boxes/Carton	Urine/Serum
HCG 116	HCG U/S Combo Cassette	GICA*	25 Tests/Box; 40 Boxes/Carton	Urine/Serum

\* GICA: Lateral Flow Gold Immunochromatographic Assay;

\*\* Carton Outer Size: 52\*45\*36CM. Gross weight about 14-15.5 KG/Carton.