


JPD-300P (7 inch)

FETAL MONITOR



JPD-300P

FETAL MONITOR (7 inch)

[Product introduction]

- The monitor of lightweight and space-saving
- It can meet the need of whole stages of labor (including antepartum and intrapartum monitoring for fetus and mother)
- It also can be connected to the obstetrical central station

[Options]

- Measurement of maternal Spo2,MHR,NIBP
- Maternal temperature measurement
- Twins monitor interface,Twins FHR measurement
- Build-in network/PC connecting interface
- FAS(fetal acoustic stimulator)

[Features]

- Lightweight, space-saving, easy operation
- Interface of parameters of FHR ,TOCO,FM,with digital and curve display
- Battery indication, pregnant No. display
- Alarm system
- 7-inch TFT color screen Folding up to 90 degree
- Build-in long life thermal printer
- DSP technology, realtime identify and measurement for FHR,accurate and reliable results
- Special high sensitive transducer
- 12 hours data storage which can be playback and printed
- Special patient event and clinical event marking bottom
- Low ultrasound power,safe to fetus
- Wide range of applying voltage(100-240V)
- Automatic Fetal Movement Detection

