

JPD-600P (12.1 inch) FETAL MATERNAL MONITOR



JPD-600P

FETAL MATERNAL MONITOR

(12.1 inch)

[Product introduction]

- The monitor of lightweight and space-saving.
- It can meet the need of whole stages of labor(including antepartum and intrapartum monitoring for fetus and mother).
- It also can be connected to the obstetrical central station



[Options]

- Twins monitor interface, Twins FHR measurement.
- Build-in network/PC/connecting interface
- Fas(fetal acoustic stimulator)



[Features]

- Light weight, compact design, space-saving, easy operation
- Interface of parameters of FHR, TOCO, FM, with digital and curve display. Maternal SpO2, HR, NIBP and temperature Maternal ECG and PR measurement.
- The setup of the system can be easy done and automatically stored.
- Alarm/system.
- 12-inch TFT color Folding up to 90 degree.
- Build-in 152 mm long life thermal printer.
- DSP technology, real time identify and measurement for FHR, accurate and reliable results.
- Special high sensitive water proof probe
- More than 24 hours data storage which can be playback and printed
- Special patient event and clinical event marking button.
- Low ultrasound power, to fetus safety
- Wide range of applying voltage (100-240V)
- Automatic Fetal Movement Detection.