

2

LEEP Surgery System



I LP−1

- LEEP Technic for uterus lesion is selected, which introduce accuracy, safety to the surgery.
- It is intended to be used for various lesion (cutting, coagulation, erosion, LEEP, Laparoscopic surgery with a variety of and gynecology surgery).
- Quick, Simple, Safety, Low Cost.
- Short Operation Time, Less Bleeding, Less Discomfort.
- Cutting wound tissue, do not affect cell morphology, which can be used for pathology diagnosis.
- Electrode can be autoclave sterilized, to ensure cutting strength and sharpness, low energy required, quick recuperate.

Model surgical mode	LP-1D	LP-1E
Pure	350W(500Ω)	350W(500Ω)
Blend 1	250W(500Ω)	250W(500Ω)
Blend 2	200W(500Ω)	200W(500Ω)
Coag 1	200W(500Ω)	200W(500Ω)
Bipo 1	70W(150Ω)	<u> </u>

