

# Vertical Pressure Steam Sterilizer

# LED display automation



#### Use

The vertical pressure steam sterilizers are assembled with heating system, micro-computer controlled system, over heat and over pressure protection system, which are reliable for the sterilizing effects, convenient of operation and energy conservation. They are ideal equipments for clinics, scientific research institutes and other organizations to sterilize surgical instruments, fabrics, glasses, culture media etc.

### Characteristics

- ◆ Fully stainless steel structure.
- ◆ LED display of working status, touch type key.
- Auto discharge the cool air, and steam discharging automatically after sterilization.
- ◆ Automatically shut off with beep reminding after sterilization.
- ◆ Easy to operate, safe and reliable.
- ◆ With two stainless steel sterilizing baskets.

## Technical Data

LS-35LD	LS-50LD	LS-75LD	LS-100LD	LS-120LD	LS-150LD
35L ф 318 × 450mm	50L ф 340 × 550mm	75L \$\phi\$ 400 \times 600mm	100L ф 440 × 650mm	120L ф 480 × 660mm	150L ф 510 × 740mm
0.22MPa					
134℃					
0.23MPa					
≤ ±1°C					
0~99min or 0~99hour59min					
105~134℃					
1×2.5KW/AC220V.50Hz	2×1.5KW/AC220V.50Hz	3×1.5KW/AC220V.50Hz		3 × 2.0KW/AC220V.50Hz	
470 × 450 × 840	520 × 500 × 970	550 × 530 × 990	580 × 560 × 1060	630 × 610 × 1130	660 × 640 × 1130
540 × 540 × 970	580 × 580 × 1090	610×610×1110	640×640×1170	690 × 690 × 1240	730 × 730 × 1230
	35L \$\dagger\$ 318 \times 450mm  \[ \begin{array}{cccccccccccccccccccccccccccccccccccc	35L \$\phi\$ 318 \times 450mm  \text{\$\frac{50L}{4340 \times 550mm}}  \text{\$\frac{1}{2}\$} \\ \begin{array}{c} \text{\$1.5KW/AC220V.50Hz} \\ \text{\$2.5KW/AC220V.50Hz} \\ \text{\$470 \times 450 \times 840}  \text{\$520 \times 500 \times 970} \end{array}	$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$	$\begin{array}{c ccccccccccccccccccccccccccccccccccc$	$\begin{array}{c ccccccccccccccccccccccccccccccccccc$