

Orbital Shaker

MPS-10



Description:

- ◆ MPS-10 Orbital Shaker is a high speed shaker with direct brushless DC motor and microprocessor controlled technique. It can be applicable in different laboratories: in microbiology, biotechnology, medical analysis, etc.

Features:

- ◆ LCD displays time and shaking speed.
- ◆ Easy to operate with one-touch knob.
- ◆ Warm root and variety of platforms can be chosen.
- ◆ When time control ends, instrument will make alarm voice.
- ◆ Gentle and reliable shaking with direct brushless DC motor.
- ◆ 12-month warranty.

Parameters:

Project	Parameters
Timing range	1min ~ 99h59min
Orbit	10mm
Capacity(PP-2 tray)	6×100ml or 4×250ml conical flask
Power	45W
Dimension(mm)	280×270×110

Project	Parameters
Shaking speed	50~300rpm
Capacity(PP-1 tray)	8pcs petri dishes
Capacity(PP-3 tray)	12×100ml or 6×250ml conical flask
Voltage	AC100V~240V, 50/60Hz
Weight(kg)	7.0