## Shaker for Microplates MS100



## **Description:**

MS100 Shaker for Microplates is a type of high speed shaker, which is with brushless DC motor and microprocessor controlled technique, specially for elisa plates(96/384 wells), tissue culture plates(24/48/96 wells).

## **Features:**

- LCD display system status and parameters.
- Stable and reliable shaking.
- Easy to operate with one touch knob.
- Automatic fault detection and buzzer alarm function.
- 12-month warranty.

## Parameters:

Project	Parameters	Project	Parameters
Timing range	1min ~ 99h59min	Shaking speed	200~1500rpm
Mixing orbit	3mm	Speed step value	10rpm
Sample capacity	4plates	Power	45W
Voltage	AC100V~240V, 50/60Hz	Dimension(mm)	280×270×110
Weight(kg)	4.0		