

# PLXF150

Mobile Surgical Table

## Application

- This operation bed for hospital treatment and surgical. Bed surface using special medical beaverboard, absorb X-ray small dose and uniform and perspective images, especially suitable for losses small C arm under X-ray TV system of inspecting and surgery

## Features

- Laminated table
- Bed surface can be take-down, moving conveniently with wheel brake control

## Main Technical Parameters

Item	Contents
Size of The Table Surface	1960mm(length)×586(width)
Table Height	660mm(lower position)~910mm(upper position)
Ascending and Descending Range	≥250mm
Load Capacity	≥135Kg

