

PLXF152

Mobile Bed for Mobile X-ray Machines

Application

- The unit is a device for the hospital's consultation, treatment and surgical use. It applies to X-ray inspection

Features

- The surface of the bed is made of medical fiberboard; The bed can rotate in 360° flexibly and the foot switch brakes quickly; The bed is low so it can satisfy the need of seriously-ill patient for special check; The longitudinal center of bed is in accord with grid box center, which ensures the symmetry of photography and filming

Main Technical Parameters

Item	Contents
Size of the Table Surface	2000mm(length)×700(width)
Table Height	500mm
Moving Range of Cassette Hanger	1680mm

