



PS T Series Ultrasonic Cleaner

Model	Vol.	Ultrasonic Power	Freq.	Time Range	Heating Temp	Heating Power	Tank Dimension	Machine Dimension	Drainage System
	L	W	KHZ	min	°C	W	(Lx w x H)mm	(Lx w x H)mm	
PS-06T	0.8	60	40	1~30	/	/	150x85x65	180x110x140	None
PS-08T	1.3	60	40	1~30	/	/	150x135x65	180x170x140	None
PS-10T	2	80	40	1~30	/	/	150x135x100	180x170x180	None
PS-20T	3	120	40	1~30	/	/	240x135x100	270x170x190	None
PS-25T	4.5	120	40	1~30	/	/	300x150x100	330x180x190	None
PS-30T	6	180	40	1~30	/	/	300x150x150	330x180x270	Yes
PS-40T	10	240	40	1~30	/	/	300x240x150	330x270x270	Yes
PS-50T	14	300	40	1~30	/	/	300x240x200	330x270x320	Yes
PS-60T	15	360	40	1~30	/	/	330x300x150	360x330x270	Yes
PS-70T	19	420	40	1~30	/	/	330x300x200	360x330x320	Yes
PS-80T	22	480	40	1~30	/	/	500x300x150	550x330x270	Yes
PS-100T	30	600	40	1~30	/	/	500x300x200	550x330x320	Yes