

Brief Introduction:

- Low pressure electro pump, power off without feedback;
- Be equipped with the liquid bottle;
- S-1 and S-2 low pressure aspirator: Alarm will occur when the fluid is full of bottle;
- S-3 has plastic mould with streamlined body; Pressure can be adjusted,1000ml liquid bottle; Alarm function will occur when the fluid is full of bottle.

Endless l ve keeps the new life safe









Standard Configuration:

Main body(including the electro pump, vacuum meter, filter, Pressure adjusting valve), liquid tube, liquid bottle.

Specification:

■ Power supply: AC220V-230V/50Hz or AC110-120V/50-60Hz

Power input: ≤ 15VAPressure range: 0-22Kpa

When the negative pressure is 22Kpa, aspirate speed: ≥ 700ml/min

■ Working noise: ≤ 55dB(A)

Capacity of liquid bottle: 1000ml (S-1)

1000ml (S-2) 1000ml (S-3)

Product Packaging:

Packing: each unit is packed into one carton

Dimension: 38.5cm x 22.5cm x 28cm

Gross weight: 4.5kg