



Pathological Sampling Bench



Configuration

UV light,
Fluorescent light,
Adjustable spotlight,
Stainless steel hanger,
"UV, FAN, LIGHT, PULVERIZER" Control.



Suction Inlet

Suction inlet adopts louver-type.
The upper mouting structure
ensures above blower suction,
avoiding air back-streaming.



Water Tap and Sink

Stainless steel
Whole-structure



Auto Flushing Design

Automatic flushing device, slightly
right-side inclined worktable, easy
for waste discharge and cleaning.



Pulverizer (Optional)

Completely smash the waste
during test, to avoid blocking
drainage pipes.

APPLICATION

Pathological sampling bench is widely used in hospital pathology department, pathology laboratory, etc. Reasonable ventilation system protects operator away from the harmful gas produced by the formalin during sampling. The hot & cold water system ensures it can adapt work in different climates.

ADVANTAGE

1. Odor-remove function.
2. Single water sink, two-person operations.
3. Exhaust air system: High-power suction device, low noise.
4. The base cabinet is inward-concave type, convenient for sitting operation.
5. Automatic flushing device, slightly right-side inclined worktable, easy for waste discharge and cleaning.
6. Optional water heater is available, realizing hot & cold water adjustment.
7. Optional pulverizer can completely smash the waste during test, to avoid blocking drainage pipes.

TECHNICAL PARAMETERS

Model	QCT-1500	QCT-1800
External Size (W*D*H)	1500*750*1950mm	1800*750*1950mm
Internal Size (W*D*H)	1395*430*900mm	1695*430*900mm
Work Surface Height	900 mm	
Consumption	300W	
Power Supply	110/220V±10%; 60/50Hz	
Exhaust Aniflow Volume	500-600m ³ /h	
Material	Main Body	304 stainless steel, thickness 1.2mm
	Work Table	304 stainless steel, thickness 1.5mm
Standard Accessories	Spotlight;	
	UV lamp;	
	Fluorescent lamp;	
	Exhaust system;	
	Automatic flushing device;	
Optional Accessories	Stainless steel hanger;	
	Adjustable hot and cold water tap;	
	180mm deep stainless steel water sink, funnel with filter.	
Optional Accessories	Pulverizer, Electrical water heater, Automatic induction hand dryer, Eye washer, Camera.	