



Standard Patient Chair
Cushion's color choose
Dark Green Sky Blue Ivory



ST-E1000
Standard Equipment

Standard Patient Chair
Cushion's color choose
Dark Green Sky Blue Ivory



Well-designed working table (area for medical instruments and medicine bottle)



Control panel



Constant temperature rinsing system

STANDARD EQUIPMENT: **OPTION**

High-grade table surface	1	Imaging system:
Liquid, Powder Spray gun(unbent3 & bent1)	4	1CCD/3CCD camera
Blowing gun	1	15 inch LCD monitor
Suction gun	1	Computer (17inch LCD display/ professional endoscope software)
Laryngoscope pre-heater	1	Printer
LED illuminating Light	1	Fiber cable
Compressor	1	Built out hernia light/built out LED cold light source
Vacuum pump	1	Endoscope:
Instrument tray	6	∅ 2.7/4mm, 0° /30° otoscope
Tweezers cup	2	∅ 2.7/4mm, 0° /30° nasoscope
Cotton cup	2	∅ 8mm, 70° /90° laryngoscope
Contaminated equipment placed tank	1	LED head light
Built-in waste tank	1	Writing table
Built-in endoscope barrel cup	1	SPECIFICATIONS:
Built-in LED cold light source	1	voltage:AC230V 50Hz
Constant temperature rinsing system	1	power:1800W
Blow-off equipment system(with warning system)	1	Body size:150*57*80 (cm)
Standard Patient Chair	1	Weight:N.W/160kg
Doctor stool	1	Weight:G.W/219kg
Microcomputer control system	1	