



### STANDARD EQUIPMENT:

①	LED illuminating light	1
②	Spray gun(unbent2 & bent1)	3
③	Suction gun	1
④	Laryngoscope pre-heater	1
⑤	Instrument tray	1
⑥	Colton cup	2
⑦	Tweezers cup	4
⑧	Medical bottle	2
⑨	Compressor	1
⑩	Vacuum pump	1
⑪	Blow-off equipment system(warning system)	1
⑫	Doctor stool	1

### SPECIFICATION:

①	Voltage:230V 50Hz
②	Power :1800W
③	Body :79*52*83cm
④	Weight :62kg
⑤	Positive pressure :1.2-2.0kg/cm <sup>2</sup>
⑥	negative pressure :700mmHg
⑦	Main suction bottle:2500cc
⑧	Pre-heater power :450W
⑨	LED illuminating Light:9V 6W



Film viewer ( Option )



Standard Patient Chair (Option)  
Cushion's color choose  
Dark Green Sky Blue Ivory

### Luxury motor patient chair technical parameters:

①	Voltage:AC230V 50Hz
②	Power :280W(max.)
③	Allowable stress:250kg
④	Seat movement itinerary:485mm(min)-625mm(max)
⑤	Hand holder is fixed backrest :90-180 degree
⑥	Pillow movement itinerary:100mm(it can fold and remove)
⑦	Weight:147kg
⑧	Seat 180degree rotate
⑨	Control by foot switch

Standard Patient Chair (Option)  
Cushion's color choose  
Dark Green Sky Blue Ivory



Luxury Motor Patient Chair(Option)  
Cushion's color choose  
Dark Green Sky Blue Ivory

### Simple motor patient chair technical parameters:

①	Voltage :AC230V 50Hz
②	Power :280W(max)
③	Allowable stress:250kg
④	Seat movement itinerary :480mm(min)-630mm(max)
⑤	Backrest 90 degree
⑥	Pillow movement itinerary :100mm(it can fold and remove)
⑦	Weight:55kg
⑧	Control by foot switch

Simple patient chair(Option)  
Cushion's color choose  
Dark Green Sky Blue Ivory



### Standard motor patient chair technical parameters:

①	Voltage:AC230V 50Hz
②	Power:280W(max)
③	Allowable stress:250kg
④	Seat movement itinerary:495mm(min)-615mm(max)
⑤	Backrest:90-135degrees
⑥	Pillow movement itinerary:100mm(it can fold and remove)
⑦	Seat 180 degree rotate
⑧	Weight:67kg
⑨	Control by foot switch