## TD4A bench top low speed centrifuge



## **Features**

1. Widely used in hospital, laboratory and blood bank.

////

- 2. Brushless DC motor for model TD4A which in great torque, free maintenance, no powder pollution, quick in speed up and down.
- 3. The flexible axle driven system which drive the rotor directly, smooth in operation, low noise and small vibration.
- 4. Micro computer control system, digital display the RCF, time and speed. There are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice.
- 5. Electric lid lock, compact design, super speed and imbalance protection.
- 6. The centrifuge body is made of high-quality steel, safe and reliable.

## **TD4A Technical Parameters:**

15-7A Toolillout I didillotoio.					
Max speed	4000rpm	Max RCF	2250xg		
Max volume	6x50ml	Noise	≤45dBA		
Timer	0~99min	Net weight	18KG		
Dimension(HxDxW)	265×483×320mm	Power supply	AC 220V 50HZ 2A		
Speed accuracy	±10rpm	Package	carton box		

## **Rotors for TD4A**

Order No.	Rotor No.	Max Speed (rpm)	Max Volume (ml)	Max RCF (g)
NO31605	Angle Rotor	4000	30×7/5ml	2250
NO31604	Angle Rotor	4000	18×10ml	2250
NO31607	Angle Rotor	4000	24×10ml	2200
NO31602	Angle Rotor	4000	12×15/7/5ml	2150
NO31603	Angle Rotor	4000	12×20ml	2220
NO31606	Angle Rotor	4000	6×50ml	2100









