



Video Laryngoscope



VL3R Reusable Blade

Safety

Vacuum packed Disposable Blades.
Screen is reinforced with anti shatter Protection.
No Cross Infection with Disposable Blades.
Anti-Fog Lens, Anti-Microbial Handle.
The Disposable blades are designed with special materials to prevent breakage, the connection part between blades and console can bear pressure more than 300kg.
316 Medical stainless steel reusable blades have the advantages of high flexure strength, no break risk; By passing through waterproof grade IPX8 testing, it can use deep water immersion disinfection.

Faster

Lightweight, Portable.



VL3D Disposable Blade

More Faster Installation and Intubation.
More Faster Removal of Blades or Drop blades.

High-Definition

The 4 inch anti-reflective display with full view has the high performance visualization capabilities for every intubation.
Better View 2.0 Mega Pixel Camera.

Easier

Rechargeable Lithium battery.
One Shot Button for Recording.

Cost Effective

Low cost disposable blades.
1000 times use 316 Medical stainless reusable blades.

Technical Specifications

Display	Size: 4 inch Full view
	Monitor: LCD
	Resolution Ratio:640X480RGB
	Aspect Ratio:4:3
	Video refresh Rata:30FPS
Data output	Easy to establish and store files
	Connect to external monitor and convenient for training And presentation
Image output	Resolution Ratio 2.0M Pixels
	Field Angle : 66°
	Illuminance > 800LUX
Laryngoscope blade Camera Include	Type Rechargeable Lithium Battery
	Long working time > 300 Minutes
	Voltage 3.7V
	Capacity 3200mAh
	Battery life Cycle: More Than 300 times
Battery	Charging times:less than 8hrs
	Temperature: -10~+45°C
	Humidity: < 93%
Transport/Storage Condition	AP 500hPa to 1060hPa
	Temperature -10~+45°C
	Humidity: 30%~85%
Working Condition	AP 700hPa to 1060hPa
	Comfortable ergonomic short handle design
Handle	

