


ELECTRONIC BALANCE



Easy to read LCD display, one touch tare up to full capacity, high resolution, quick readable weighting, easy to operate, widely used in high precise industrial applications, laboratory, university.

MODEL	XY600-1B	XY1000-1B	XY2000-1B	XY3000-1BF	XY5000-1BF
Capacity (g)	610	1100	2100	3100	5100
Price(USD)					
Resolution (g)	0.1				
Min weighing (g)	0.4				
Stable Time	≤ 2 S				
Pan Size	Φ 130 mm			160X160mm	
Operation temp	5-25℃				
Repeat ability	0.2g				
Liner	0.3g				
Cal.weight	200g	500g	1000g	2000g	
Option	RS232C & RS485、 Printer				
Packing Size	Inside:310×230×130mm			Outside:660×470×320mm	
G.weight	2.0/22kg				

* The calibration weight is optional parts when the capacity ≥2000g or resolution >0.01g.