



The LED display are suit for most use place. XY-C series is widely used in industrial applications ,laboratory and university because of it made of stable quality.We can use the under weighing functions to test the density.

MODEL	XY200C	XY300C	XY500C	XY600C	XY1000C	XY2000C	XY3000C
Capacity (g)	210	310	510	610	1100	2100	3100
Price(USD)							
Resolution (g)	0.01						
Min weighing (g)	0.04						
Stable Time	≤ 2 S						
Pan Size	Φ 130 mm						
Operation temp	5-25℃						
Repeat ability	0.02g						
Liner	0.03g						
Cal.weight	200g	200g	200g	200g	500g	1000g	2000g
Option	RS232C & RS485、 Printer						
Packing Size	Inside:305×215×115mm			Outside:660×470×320mm			
G.weight	2.0/22kg						


* The calibration weight is optional parts when the capacity ≥2000g or resolution >0.01g.

PRECISION BALANCE

Standard Configuration:



Optional Accessories:



XY-C