



XY -MA/MB series weighing scale use high precision weighing sensor to keep weighing stable. we support different interface and accessories to extend functions. and the under weighing function is suit for density test.

WEIGHING SCALE

XY-MA/MB Series

MODEL	XY10MA	XY15MA	XY20MA	XY30MA	XY5MB	XY10MB	XY15MB	XY20MB	XY25MB
Capacity	11kg	16kg	21kg	31kg	5.5kg	11kg	16kg	21kg	28kg
Price(USD)									
Sensor	Load Cell								
Resolution (g)	1				0.1				
Min weighing (g)	4				0.4				
Stable Time	≤ 2 S								
Display	LCD (white back light with black font)								
Pan Size	315x230mm				315x230mm				
Operation temp	5-25°C								
Repeat ability	±2g				±0.2g				
Liner	±3g				±0.3g				
Draft shield size	240mmx190mmx265mm								
Option	RS232C & RS485、 Printer								
Packing Size	Inside:410×165×400mm				Outside:800×420×420mm				
G.weight	4/16 kg								

* The calibration weight is optional parts when the capacity ≥2000g or resolution >0.01g.

Standard Configuration:



Optional Accessories:

