


\Rightarrow

PORTABLE OXYGEN CONCENTRATOR




Features

- Compact design, light weight and easy to carry with trolley bag.
- 1-7L flow adjustable, with digital oxygen flow meter.
- Real-time monitoring oxygen concentration and flow, safer and more reliable.
- Innovative cooling technology, ultra quiet.

DIMENSION

- Powered by AC 220V/110V and vehicle DC 12V.
- For pregnant women, sub-healthy people and patient short time outdoor oxygen supplement.





PRODUCT INFORMATION

FLOW RATE	1–7L
OUTLET PRESSURE	0.04-0.07Mpa
OXYGEN CONCENTRATION	93%–38%
OUTLET PRESSUER	6-10Psi (0.04-0.07Mpa)
NET WEIGHT	6kgs

219*276*313mm