

Tunnel-type electronic blood pressure monitor

Protect your parents blood pressure health



Simplify and measure with one click
Measure directly without a cuff

Reach out and measure without assistance to avoid inaccurate measurement caused by improper wearing
It is more suitable for middle-aged and elderly people



Full voice guidance

Parents can use it easily
without any off-site assistance

||| High pressure: 120mmHg

||| Low pressure: 80mmHg

||| The pulse is 80 beats
per minute



The arm rest is foldable
It doesn't take up
much space

More support points ensure
more accurate pressure
measurement



Give you a comfortable experience



One-click operation

One-click intelligent pressurization
Achieve rapid measurement

Elbow pressure switch

The test requires the elbow area
Only by holding the switch can one succeed
Start measuring



In case of emergency
One-click emergency stop

In case of emergency One-click emergency stop

Storage of multiple sets of blood pressure data
Clearly display the trend changes



Product parameters



Product Name: Arm Electronic Blood Pressure Monitor

Product model: KWL-819

Display mode: LED

Measurement method: oscilloscope method

The main unit weighs approximately 1268 grams (excluding the battery).

Operating voltage: 6V d.c. (4 AA batteries)

The measurable arm circumference is approximately 22 to 32cm

Pressure measurement error: $\pm 3\text{mmHg}(\pm 0.4\text{kPa})$

The measurement range of pulse rate is 40 to 160 beats per minute

The main unit measures 257.8mm in length, 202mm in width and 215mm in height

The screen size is 100mm*93mm