

## **INDIRECT OPHTHALMOSCOPE**

## YZ25C

## **Indirect Ophthalmoscope (Rechargeable)**

- LED lamp(Low-temperature,Low-pressure,High-brightness),the light is uniformly,bright and brightness can adjust continuously.
- Imaging clearly,strong dimensional sense and large field of view range. The main body part of this device adopts engineering plastics, compact and exquisite.
- Super lightweight design to reduce fatigue, operation is very convenient.
- The portable power is suitable for surgery,ward round and consultation.
- Equipped with color filter, teaching mirror and sclera depressor, operating be more handy.





Specifications	
Angle of View	45° ±5°
Pupil Distance	52mm ~ 74mm
Weight	700g
Adjustable Range of the Head Frame	520mm ~ 640mm
Depth	85mm ~ 125mm
Illumination Source	LED
Power Supply	Lit battery pack DC:7.4V Rechargeable